

東ソーマルチコントロール（心疾患）参考値一覧

2014. 12. 09 updated

lot	BNP (pg/mL)		cTnI3 (ng/mL)		CK-MB (ng/mL)		MYO (ng/mL)	
	level1	level2	level1	level2	level1	level2	level1	level2
E349635	33-49	390-590	0.63-0.94	6.0-9.1	8.8-13	91-140	42-63	220-330
E649636	34-51	380-570	0.61-0.91	5.8-8.7	8.5-13	88-130	42-63	230-340
E949637	33-49	380-560	0.41-0.62	4.0-6.0	8.2-12	83-120	41-61	220-330
EY49638	34-51	330-500	0.44-0.67	4.7-7.1	7.8-12	70-110	43-64	230-350