

東ソーコントロールセット（心疾患）参考値一覧

2015.09.16 updated

Lot.	BNP (pg/mL)		cTnI3 (ng/mL)		CK-MB (ng/mL)		MYO (ng/mL)		CL-BNP (pg/mL)	
	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2
E949637	33-49	380-560	0.41-0.62	4.0-6.0	8.2-12	83-120	41-61	220-330		
EY49638	34-51	330-500	0.44-0.67	4.7-7.1	7.8-12	70-110	43-64	230-350	31-46	300-450
F349639	32-49	380-570	0.39-0.59	4.0-6.0	7.5-11	79-120	28-42	210-310	31-46	330-500
F749640	35-52	380-570	0.43-0.64	4.4-6.6	7.5-11	84-130	41-61	220-340	32-48	330-490

Lot.	使用期限
E949637	2015.09
EY49638	2015.11
F349639	2016.03
F749640	2016.07