

東ソーコントロールセット（心疾患）参考値一覧

2017.11.14 updated

Lot.	BNP (pg/mL)		cTnI3 (ng/mL)		CK-MB (ng/mL)		MYO (ng/mL)		CL-BNP (pg/mL)		CL-cTnI (pg/mL)		CL-CKMB (ng/mL)		CL-MYO (ng/mL)	
	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2
GY49644	33-49	380-570	0.36-0.54	3.8-5.7	6.8-10	81-120	38-57	220-340	31-46	330-490	250-380	2800-4200				
H249645	35-52	380-570	0.41-0.61	4.3-6.5	8.6-13	78-120	46-69	240-360	30-45	320-480	300-450	3500-5200	7.8-12	68-100	43-64	210-320
H649646	34-52	380-580	0.46-0.70	4.1-6.2	8.2-12	86-130	42-63	230-340	30-46	320-490	360-540	3300-5000	7.4-11	74-110	39-59	220-330
HX49647	34-50	380-580	0.39-0.58	4.2-6.4	9.2-14	88-130	39-58	230-350	30-45	320-480	250-370	3000-4500	8.1-12	76-110	34-50	210-310

Lot.	使用期限
GY49644	2017.11
H249645	2018.02
H649646	2018.06
HX49647	2018.10