

東ソーコントロールセット（心疾患）参考値一覧

2018.11.13 updated

| Lot. | BNP (pg/mL) | | cTnI3 (ng/mL) | | CK-MB (ng/mL) | | MYO (ng/mL) | | CL-BNP (pg/mL) | | CL-cTnI (pg/mL) | | CL-CKMB (ng/mL) | | CL-MYO (ng/mL) | |
|---------|-------------|---------|---------------|---------|---------------|--------|-------------|---------|----------------|---------|-----------------|-----------|-----------------|--------|----------------|---------|
| | level1 | level2 | level1 | level2 | level1 | level2 | level1 | level2 | level1 | level2 | level1 | level2 | level1 | level2 | level1 | level2 |
| I249648 | 31-46 | 400-610 | 0.39-0.58 | 4.3-6.4 | 7.2-11 | 73-110 | 39-58 | 220-330 | 27-40 | 330-500 | 280-420 | 3500-5200 | 6.5-9.7 | 66-99 | 35-52 | 210-310 |
| I649649 | 34-50 | 400-600 | 0.38-0.56 | 4.2-6.3 | 7.4-11 | 80-120 | 36-53 | 210-320 | 28-42 | 340-510 | 270-400 | 3100-4700 | 6.9-10 | 74-110 | 32-48 | 200-290 |
| I949650 | 32-48 | 400-600 | 0.37-0.56 | 4.0-6.0 | 7.6-11 | 78-120 | 45-67 | 220-330 | 27-40 | 330-490 | 250-370 | 2800-4200 | 6.9-10 | 70-110 | 39-58 | 200-300 |

| Lot. | 使用期限 |
|---------|---------|
| I249648 | 2019.02 |
| I649649 | 2019.06 |
| I949650 | 2019.09 |