

東ソーコントロールセット（心疾患）参考値一覧

2016.04.04 updated

Lot.	BNP (pg/mL)		cTnI3 (ng/mL)		CK-MB (ng/mL)		MYO (ng/mL)		CL-BNP (pg/mL)	
	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2
F749640	35-52	380-570	0.43-0.64	4.4-6.6	7.5-11	84-130	41-61	220-340	32-48	330-490
FY49641	35-53	390-590	0.38-0.57	4.2-6.3	8.2-12	88-130	37-56	220-340	31-47	340-510
G349642	32-49	360-550	0.38-0.58	4.1-6.1	8.1-12	87-130	29-43	210-320	31-46	330-490

Lot.	使用期限
F749640	2016.07
FY49641	2016.11
G349642	2017.03