

## 東ソーコントロールセット（心疾患）参考値一覧

2017. 3. 15 updated

Lot.	BNP (pg/mL)		cTnI3 (ng/mL)		CK-MB (ng/mL)		MYO (ng/mL)		CL-BNP (pg/mL)		CL-cTnI (pg/mL)		CL-MYO (ng/mL)		CL-CKMB (ng/mL)	
	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2
G349642	32-49	360-550	0.38-0.58	4.1-6.1	8.1-12	87-130	29-43	210-320								
G649643	33-50	370-560	0.36-0.54	4.0-5.9	7.6-11	77-120	39-58	220-330	30-44	320-480	240-360	2700-4100				
GY49644	33-49	380-570	0.36-0.54	3.8-5.7	6.8-10	81-120	38-57	220-340	31-46	330-490	250-380	2800-4200				
H249645	35-52	380-570	0.41-0.61	4.3-6.5	8.6-13	78-120	46-69	240-360	30-45	320-480	300-450	3500-5200	43-64	210-320	7.8-12	68-100

Lot.	使用期限
G349642	2017. 03
G649643	2017. 06
GY49644	2017. 11
H249645	2018. 02