

東ソーコントロールセット（心疾患）参考値一覧

2020.7.7 updated

Lot.	BNP (pg/mL)		cTnI3 (ng/mL)		CK-MB (ng/mL)		MYO (ng/mL)		CL-BNP (pg/mL)		CL-cTnI (pg/mL)		CL-CKMB (ng/mL)		CL-MYO (ng/mL)	
	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2
JX49653	29-43	390-580	0.37-0.56	3.9-5.9	7.9-12	87-130	41-62	220-330	27-41	340-510	260-390	3100-4600	7.3-11	77-110	36-55	200-300
A249654	31-47	310-470	0.33-0.50	3.5-5.3	7.1-11	75-110	39-59	190-290	28-43	270-410	220-330	2500-3800	6.7-10	70-100	35-53	180-270
A649655	30-44	380-580	0.40-0.60	4.3-6.5	7.4-11	82-120	40-60	230-340	27-41	340-510	300-450	3300-4900	7.4-11	75-110	37-56	210-310

Lot.	使用期限
JX49653	2020.10
A249654	2021.02
A649655	2021.06