

東ソーコントロールセット（心疾患）参考値一覧

2016.10.17 updated

| Lot. | BNP (pg/mL) | | cTnI3 (ng/mL) | | CK-MB (ng/mL) | | MYO (ng/mL) | | CL-BNP (pg/mL) | | CL-cTnI (pg/mL) | |
|---------|-------------|---------|---------------|---------|---------------|--------|-------------|---------|----------------|---------|-----------------|-----------|
| | level1 | level2 | level1 | level2 | level1 | level2 | level1 | level2 | level1 | level2 | level1 | level2 |
| FY49641 | 35-53 | 390-590 | 0.38-0.57 | 4.2-6.3 | 8.2-12 | 88-130 | 37-56 | 220-340 | | | | |
| G349642 | 32-49 | 360-550 | 0.38-0.58 | 4.1-6.1 | 8.1-12 | 87-130 | 29-43 | 210-320 | | | | |
| G649643 | 33-50 | 370-560 | 0.36-0.54 | 4.0-5.9 | 7.6-11 | 77-120 | 39-58 | 220-330 | 30-44 | 320-480 | 240-360 | 2700-4100 |

| Lot. | 使用期限 |
|---------|---------|
| FY49641 | 2016.11 |
| G349642 | 2017.03 |
| G649643 | 2017.06 |