

東ソーコントロールセット（心疾患）参考値一覧

2021.03.08 updated

Lot.	BNP (pg/mL)		cTnI3 (ng/mL)		CK-MB (ng/mL)		MYO (ng/mL)		CL-BNP (pg/mL)		CL-cTnI (pg/mL)		CL-CKMB (ng/mL)		CL-MYO (ng/mL)	
	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2
A649655	30-44	380-580	0.40-0.60	4.3-6.5	7.4-11	82-120	40-60	230-340	27-41	340-510	300-450	3300-4900	7.4-11	75-110	37-56	210-310
A949656	33-49	400-590	0.41-0.61	4.6-7.0	7.4-11	86-130	41-61	240-360	29-44	340-510	280-410	3400-5200	7.0-11	78-120	33-50	210-310
B149657	32-49	390-590	0.41-0.62	4.7-7.1	8.3-12	92-140	43-64	240-360	28-43	340-510	290-440	3600-5500	7.5-11	82-120	41-62	230-340

Lot.	使用期限
A649655	2021.06
A949656	2021.09
B149657	2022.01