

東ソーマルチコントロール（心疾患）参考値一覧

2013.02.12update

lot	BNP(pg/ml)		cTnI3(ng/ml)		CK-MB(ng/ml)		MYO(ng/ml)	
	level1	level2	level1	level2	level1	level2	level1	level2
B749621	30-45	400-600	0.62-0.93	5.8-8.7	6.5-9.7	69-100	29-43	220-330
B949622	32-49	370-560	0.63-0.94	5.8-8.7	7.3-11	75-110	45-67	230-340
BY49623	32-49	370-560	0.60-0.90	5.7-8.5	7.2-11	68-100	37-55	230-350
C149624	34-51	380-570	0.65-0.97	5.7-8.6	7.5-11	78-120	46-70	250-370
C349625	32-48	380-570	0.55-0.82	5.5-8.3	7.4-11	73-110	41-61	200-310
C549626	34-51	370-550	0.58-0.87	5.6-8.4	6.7-10	75-110	42-63	230-350
C749627	32-48	370-560	0.63-0.94	5.9-8.8	7.3-11	77-120	40-59	220-330
C949628	32-48	390-580	0.58-0.87	5.2-7.7	6.7-10	70-100	39-58	240-370
CY49629	32-49	400-600	0.59-0.88	6.3-9.4	7.0-10	71-110	40-59	240-360