

## 東ソーコントロールセット（心疾患）参考値一覧

2019.9.11 updated

Lot.	BNP (pg/mL)		cTnI3 (ng/mL)		CK-MB (ng/mL)		MYO (ng/mL)		CL-BNP (pg/mL)		CL-cTnI (pg/mL)		CL-CKMB (ng/mL)		CL-MYO (ng/mL)	
	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2	level1	level2
I949650	32-48	400-600	0.37-0.56	4.0-6.0	7.6-11	78-120	45-67	220-330	27-40	330-490	250-370	2800-4200	6.9-10	70-110	39-58	200-300
J249651	29-44	400-590	0.35-0.52	4.0-6.0	7.3-11	79-120	44-66	230-360	25-38	330-490	260-390	3200-4800	6.8-10	71-110	37-55	200-300
J649652	30-45	380-570	0.36-0.54	3.9-5.9	7.7-12	82-120	41-61	230-350	26-39	330-490	240-350	2900-4400	7.1-11	74-110	34-51	210-310

Lot.	使用期限
I949650	2019.09
J249651	2020.02
J649652	2020.06